

Download Hyfrecator 2000 Service Manual

Page 1. Operator's Manual 2000 Serving the Physician Since 1937... Page 2. For Technical Services or Return Authorization Phone: 303-699-7600 / 1-800-552-0138 Extension 5274 Fax 303-699-9854 For Customer Services or to order parts phone: 1-800-448-6506 / 315-797-8375 / Fax 315-735-6235 or contact your CONMED Representative. The Hyfrecator ® is a multi-function electrosurgical generator with both monopolar and bipolar modes, self-diagnostics, and recall of stored power settings. Upon installation physicians experience simple use, precise tissue therapy, and years of safe, reliable operation. There is only one Hyfrecator®. The Hyfrecator ® is a multi-function electrosurgical generator with both monopolar and bipolar modes, self-diagnostics, and recall of stored power settings. Upon installation physicians experience simple use, precise tissue therapy, and years of safe, reliable operation. There is only one Hyfrecator®. MFR#: CMD-7-900-OM-ENG. Conmed Hyfrecator 2000 Operator's Manual in English, 7-900-OM-ENG. Includes safety precautions, terminology for systems and procedural configurations, applications, maintenance, optional accessories, specifications and performance.